



# PCMS01A-HC



## Product

Patient Care Monitoring System

### Description:

Four-Digit Display with 1" High Digits and Red LEDs, Four-Digit Display with 1" High Digits, Red-Yellow-Green lamps with 2.5" x 2.5" LED Clusters, in Common Metal Enclosure for Hospital Compliance.

**Item:** PCMS01A-HC

**UPC:** 788581652824

**EAN:** 0788581652824

**Stock Status:** Special order item. Allow 2 weeks lead time.

**Dimensions:** 21.5" x 8" x 2.5"

**Weight:** 6 lbs.

**MSRP:** \$995.00

## Specifications

- Two RJ45 modular connectors for RS422 serial data transmission and reception.
- RJ11 modular connector for remote Start, Stop and Reset inputs.
- Power connector (2.1x5.5mm).
- Power requirements: 10-15VAC or 12-20VDC.
- Current consumption: 500mA maximum.
- Construction: Stainless steel metal enclosure with overlay having transparent red, yellow and green windows.

### Mounting methods

- Standard mounting fits over a single or double gang electrical box for wall mounting.
- Alternate mounting is 2 key hole slots for wall mounting.

## Features

- Four digit LED display showing the number of minutes and seconds since the patient was last checked on.
- Four digit LED display showing the number of hours and minutes that the patient has been in the room.
- Green indicator indicate that the patient has recently been checked on.
- Yellow indicator indicate that the patient needs to be checked on soon.
- Red indicator to indicate that the patient is over due to be checked on.
- User programmable times for the Green, Yellow and Red indicator lamps using an IR remote.

## In the Box

- PCMS01A Dual LED display with RYG Indicators
- IR27A Infrared remote.
- RC Reference guide
- WD Wiring diagram
- M8BK003F Cable, 3ft, Cat-5
- WT12VAC1200MA-RA 12V AC, 1200mA transformer
- KP2-NTGR Plate with no-touch button and reset button